



Manufacturing Innovative IVD for the World Since 1973

# ELISA & LIA

KITS



ELISA KIT NAME	REF	ASSAY TIME	SENSITIVITY	STABILITY
5 $\alpha$ -Androstane 3 $\alpha$ , 17 $\beta$ -Diol Glucuronide (3 $\alpha$ -Diol G)	CAN-DG-460	45 minutes	0.1 ng/mL	18 months
25-Hydroxyvitamin D [25(OH)D]	CAN-VD-510	105 minutes	5.5 ng/mL	18 months
Adiponectin	CAN-APN-5000	105 minutes	0.06 ng/mL	18 months
Aldosterone	CAN-ALD-500	80 minutes	9.1 pg/mL	18 months
Androstenedione	CAN-AD-208	80 minutes	0.04 ng/mL	18 months
Anti-Thyroglobulin (Tg)	CAN-ATG-4150	80 minutes	2.06 IU/mL	12 months
Anti-Thyroid Peroxidase (TPO)	CAN-ATP-4160	80 minutes	0.819 IU/mL	12 months
$\beta$ 2-Microglobulin	CAN-B-4300	80 minutes	0.1 mg/L	18 months
C-Peptide	CAN-CP-4380	110 minutes	0.2 ng/mL	18 months
C-Reactive Protein, High Sensitivity (hs-CRP)	CAN-CRP-4360	60 minutes	10 ng/mL	18 months
Chorionic Gonadotropin (hCG)	CAN-HCG-4120	75 minutes	0.7 IU/L	12 months
Cortisol Saliva	CAN-C-240	65 minutes	0.033 ng/mL	12 months
Cortisol Saliva	CAN-C-290	65 minutes	1 ng/mL	12 months
Cortisol Serum	CAN-C-270	60 minutes	2.061 $\mu$ g/dL	18 months
Dehydroepiandrosterone (DHEA)	CAN-DH-490	110 minutes	0.092 ng/mL	18 months
Dehydroepiandrosterone Sulfate (DHEAS)	CAN-DHS-480	65 minutes	0.005 $\mu$ g/mL	18 months
Dihydrotestosterone (DHT)	CAN-DHT-280	75 minutes	6 pg/mL	12 months
Estradiol	CAN-E-430	75 minutes	10 pg/mL	18 months
Estriol, Free	CAN-E-640	80 minutes	0.058 ng/mL	12 months
Estriol, Free – Serum & Saliva	CAN-E-620	80 minutes	0.058 ng/mL 0.034 ng/mL	12 months
Estrogens, Total	CAN-E-630	180 minutes	12.4 pg/mL	12 months
Estrone	CAN-E-420	80 minutes	14.8 pg/mL	18 months
Ferritin	CAN-F-4280	45 minutes	0.44 ng/mL	18 months
Follicle Stimulating Hormone (hFSH)	CAN-FSH-4060	80 minutes	1 IU/L	18 months
Growth Hormone (hGH)	CAN-GH-4070	75 minutes	0.2 ng/mL	18 months
Insulin-Like Growth Factor Binding Protein-1 (IGFBP-1)	CAN-IGF-4140	75 minutes	0.5 $\mu$ g/L	18 months
Leptin	CAN-L-4260	105 minutes	0.5 ng/mL	18 months
Luteinizing Hormone (hLH)	CAN-LH-4040	80 minutes	0.2 IU/L	18 months
Plasma Renin Activity (PRA)	CAN-RA-4600	105 minutes	0.166 ng/mL of Angiotensin-I	12 months
Pregnenolone	CAN-PRE-4500	75 minutes	0.05 ng/mL	18 months
Progesterone Saliva	CAN-P-310	80 minutes	20 pg/mL	12 months
Progesterone Serum	CAN-P-305	80 minutes	0.1 ng/mL	18 months
Progesterone, 17 $\alpha$ -OH	CAN-P-400	80 minutes	0.051 ng/mL	18 months

ELISA KIT NAME	REF	ASSAY TIME	SENSITIVITY	STABILITY
Prolactin	CAN-PRL-4100	75 minutes	10 µIU/mL	18 months
Prostate Specific Antigen, Free (fPSA)	CAN-FPSA-4400	75 minutes	0.05 ng/mL	18 months
Prostate Specific Antigen, Total (tPSA)	CAN-TPSA-4300	105 minutes	0.1 ng/mL	12 months
Resistin	CAN-RSN-4000	105 minutes	0.004 ng/mL	12 months
Sex Hormone Binding Globulin (SHBG)	CAN-SHBG-4010	60 minutes	0.1 nmol/L	18 months
Testosterone, Free	CAN-FTE-260	75 minutes	0.018 pg/mL	18 months
Testosterone, Total	CAN-TE-250	75 minutes	0.022 ng/mL	18 months
Testosterone Saliva	CAN-TE-300	90 minutes	6.75 pg/mL	12 months
Thyroid Stimulating Hormone (TSH)	CAN-TSH-4080	105 minutes	0.1 µIU/mL	18 months
Thyroxine, Free (fT4)	CAN-FT4-4340	75 minutes	1 pg/mL	18 months
Thyroxine, Total (T4)	CAN-T4-4240	50 minutes	0.6 µg%	18 months
Triiodothyronine, Free (fT3)	CAN-FT3-4230	75 minutes	0.3 pg/mL	18 months
Triiodothyronine, Reverse (rT3)	CAN-RT3-100	105 minutes	0.009 ng/mL	12 months
Triiodothyronine, Total (T3)	CAN-T3-4220	75 minutes	0.16 ng/mL	18 months

LIA KIT NAME	REF	ASSAY TIME	SENSITIVITY	STABILITY
5α-Androstane 3α, 17β-Diol Glucuronide (3α-Diol G)	CAN-DG-6010	60 minutes	0.1 ng/mL	12 months
Cortisol	CAN-C-7010	90 minutes	0.15 µg/dL	12 months
Dehydroepiandrosterone Sulphate (DHEAS)	CAN-DHS-7020	60 minutes	0.2 µg/mL	12 months
Estriol Saliva	CAN-E-7030	95 minutes	0.03 ng/mL	12 months
Estrone	CAN-E-6060	90 minutes	8.8 pg/mL	12 months
Estrone Saliva	CAN-E-7000	90 minutes	1.2 pg/mL	12 months
Plasma Renin Activity (PRA)	CAN-RA-7070	120 minutes	0.22 ng/mL of Angiotensin-I	12 months
Progesterone	CAN-P-6050	90 minutes	0.1 ng/mL	12 months
Prostate Specific Antigen	CAN-TPSA-6020	120 minutes	0.082 ng/mL	12 months
Testosterone, Free	CAN-FTE-6080	90 minutes	0.17 pg/mL	12 months
Testosterone Saliva	CAN-TE-6090	120 minutes	1 pg/mL	12 months
Thyroid Stimulating Hormone	CAN-TSH-6030	90 minutes	0.13 µIU/mL	12 months
Thyroxine, Free (fT4)	CAN-FT4-6040	90 minutes	1 pg/mL	12 months



For more information, please visit

**dbc-labs.com**

